

Universal Frame Outer Socket for Trans-Femoral Amputees (UFOS™-TF)

The UFOS™-TF, designed to be used in conjunction with either the FLO-TECH-TOR™-TF, permits controlled, measurable weight bearing, as part of the APOPPS®-TF. This combination enables the patient to gradually build a tolerance to the firm support required in a properly fitted prosthesis. During this critical period of rehabilitation..

As the patient progresses, the semi-rigid plastic sleeve promotes good venous return and helps shape and prepare the residual limb for a permanent prosthesis, while helping restore normal gait and stance.

The UFOS™-TF is a one-piece preparatory outer socket that allows the Prosthetist to attach the components of choice. When used with either the FLO-TECH-TOR™-TF offers opportunity for ambulatory training from post-operative amputation to definitive prosthesis.

The UFOS™-TF with foot and pylon components can be removed without removing the inner socket, offering practitioners a versatile rehabilitation system. The patient can be sent home, or back to their hospital room with the FLO-TECH-TOR™-TF socket on. The UFOS™-TF remains in the control and possession of the Physical Therapist or Prosthetist until the patient is ready to assume the responsibility of donning and doffing it themselves.

Each UFOS™-TF is sized to fit with the corresponding FLO-TECH-TOR™-TF your patient is using. Sizes will fit ischial circumference measurements of 16" to 32", standard sockets are 12" in length. We offer Prosthetists 4 pre-fabricated preparatory trans-femoral UFOS™-TF sockets available, to fit many trans-femoral amputees. Sockets include all straps, bands and foam distal end pads.

Custom sockets can also be fabricated from your measurements. We make every reasonable effort to ship Custom orders within 24-48 hours, making the benefits of the UFOS™-TF available to most of trans-femoral amputees.



FLO-TECH® O & P Systems, Inc.
PO Box 462 - 7325 Halseyville Road
Trumansburg, New York 14886
www.1800FLO-TECH.com

All materials © FLO-TECH® October 2017

Phone: 1-800-FLO-TECH (356-8324)
Phone: (607) 387-3070
Fax: (607) 387-3176
E-mail: info@1800FLO-TECH.com