

RESIDUUM: LEFT ____ RIGHT ____ 29 DEGREES OR LESS ____ 30 DEGREES OR MORE ____

MEASUREMENTS (ALL MEASUREMENTS TAKEN OVER BANDAGES, SOCKS, ETC.):

LENGTH: MPT TO DISTAL END ____ In. (Cm) (shorter than 4" or longer than 8 1/4" without bandages requires custom made sockets)

CIRCUMFERENCES: MPT ____ In. (Cm) 2" UP FROM DISTAL END ____ In. (Cm) (Choose size based on this Circ.)

Sockets are automatically sized to provide the ideal anatomic knee circumference for the socket size ordered.

ADULT SIZES

#1416 – (MPT 14" TO 16")

Distal: 13" (33cm) to 15" (38cm)

Range: 11" (28cm) to 16" (41cm)

MODIFIED FLO-TECH-TOR™

L ____ 7 or 9

R ____

29 DEGREES OR LESS:

The anterior distal portion of a FLO-TECH-TOR™ follows the angle of the contractures (up to 29 degrees) while the posterior section allows for graduated reduction of the contracture.

#1618 – (MPT 16" TO 18")

Distal: 15" (38cm) to 17" (43cm)

Range: 13" (33cm) to 18" (46cm)

L ____ 7 or 9

R ____

#1821 – (MPT 18" TO 21")

Distal: 17" (43cm) to 19" (48.5cm)

Range: 15" (38cm) to 22" (56cm)

L ____ 7 or 9

R ____

*30 DEGREES OR MORE:

A modified thigh section of a FLO-TECH-TOR™ & complete VCSPS™ with jointed uprights.

#2124 – (MPT 21" TO 24")

Distal: 20" (51cm) to 23" (58.5cm)

Range: 17" (43cm) to 25" (63.5cm)

L ____ 7 or 9

R ____


These socket systems include:

- An additional neoprene band with Velcro® type attachment at the knee
- A long nylon stockinet pull-sock to allow for periodic re-positioning of tissue
- Slots cut in the distal pads
- Distal hole in the socket

*Create a jointed APOPPS® (Rehab and Prep prosthesis all in one) by adding an UFOS™ and waist belt, which allows for the attachment of Pylon and Foot components for ambulatory training. Remove the joints and thigh section (from the flexion control socket) when the contracture is eliminated to create a Prep System.

COMFORT OPTIONS & COMPONENT ITEM PRICES:

QTY	DESCRIPTION	QTY	DESCRIPTION
_____	Modifications to standard sockets	_____	Each Amputation Sock SIZE: _____
_____	Waist belt/suspension strap	_____	Post-Op Foot - L <input type="checkbox"/> R <input type="checkbox"/> Size: _____
_____	Fork/suspension strap	_____	Titanium SACH Foot Adapter
_____	Per Socket Strap - ALL SIZES	_____	Titanium Pyramid Socket Attachment Plate
_____	Velcro® brand Compatible Neoprene socket band _____	_____	Titanium Pylon Tube w/adaptor
_____	Distal end pads - Per Pkg. SIZE: _____	_____	Titanium Tube Clamp
_____	Distal end pad, Reticulated SIZE: _____	_____	S.S. SACH Foot Adapter
_____	PADDED/LINED Waist belt	_____	S.S. Pyramid Socket Attachment Plate
_____	Pour & Modify Cast	_____	S.S. Pylon Tube w/adaptor
_____	4 bolt adapter plate, 4 short, 4 long screws	_____	S.S. Tube Clamp



FLO-TECH® O & P Systems, Inc.
 PO Box 462 - 7325 Halseyville Road
 Trumansburg, New York 14886
 www.1800FLO-TECH.com

All materials © FLO-TECH® Modified March 2017

Phone: 1-800-FLO-TECH (356-8324)

Phone: (607) 387-3070

Fax: (607) 387-3176

E-mail: info@1800FLO-TECH.com