

# Trans-Femoral APOPPS® -TF

# Measurement/Order Form

PROSTHETIST: \_\_\_\_\_ PO NUMBER: \_\_\_\_\_ Date: \_\_\_\_\_

Use the following sizing information for selecting APOPPS®-TF sockets for your patient. Each of these sockets is a one-piece, polyethylene, clamshell design with anatomically trimmed overlapping sides, ideal for accommodating bulbous, straight or conical shapes and gently reducing the residuum throughout the rehabilitation period.

- ⊗ The FLO-TECH-TOR™-TF is used as a prosthetic socket and for early supervised weight bearing prosthetic management of trans-femoral amputations.
- ⊗ The UFOS™-TF is an external frame outer socket used to attach the pylon and foot unit. It is sized to match the FLO-TECH-TOR™-TF you select for your patient. The APOPPS®-TF Program allows your patient to wear the FLO-TECH-TOR™ socket alone when the pylon and foot are not needed.

**The APOPPS® Program is designed to maintain a healthy amputation throughout the rehabilitation process.**

**MEASUREMENTS** (ALL MEASUREMENTS SHOULD BE TAKEN OVER BANDAGES, SOCKS, ETC.): **RESIDUUM: Left** \_\_\_\_\_ **Right** \_\_\_\_\_

**LENGTH Ischium to Distal End:** \_\_\_\_\_ in/cm **CUSTOM** \_\_\_\_\_

**\*CIRCUMFERENCES:** Ischial Level \_\_\_\_\_ in/cm **\*\*3" Up from Distal End** \_\_\_\_\_ in/cm

*\*Top measurement at each size indicates ideal Ischial level circumference*  
*\*\*Bottom measurement indicates the minimum to maximum circumference range at the Distal level.*  
*Sockets are automatically sized to fit the ideal circumference at the Above Ischial location.*

12" (9¾" to 11¼" LENGTH)	FLO-TECH-TOR™	UFOS™	APOPPS®
<b>SMALL – *16" (41cm) to 22" (56cm)</b> **Range: 13" (33cm) to 17" (43cm)	L _____ R _____	L _____ R _____	L _____ R _____
<b>MEDIUM – *19" (48cm) to 25" (63.5)</b> **Range: 16" (41cm) to 20" (51cm)	L _____ R _____	L _____ R _____	L _____ R _____
<b>LARGE – *22" (56cm) to 28" (71cm)</b> **Range: 19" (48cm) to 23" (58.5cm)	L _____ R _____	L _____ R _____	L _____ R _____
<b>X-LARGE – *26" (66cm) to 32" (81.5cm)</b> **Range: 24" (61cm) to 28" (71.5cm)	L _____ R _____	L _____ R _____	L _____ R _____

**Knee, Foot & Pylon Units Available:** Includes manual locking knee, rotatable 4-hole connector, S.A.C.H. foot adapter, pylon and **Post-Op Wayfarer Foot** L  R  Size: 24  26  28  (Required with all foot orders)  
 Titanium/Knee/Pylon Unit  Stainless Steel/Knee/Pylon Unit

**COMFORT OPTIONS & COMPONENT ITEMS:**

QTY	DESCRIPTION	QTY	DESCRIPTION
_____	Modification to Standard Socket	_____	Each Amputation Socks Size: _____
_____	Neoplush Pelvic Suspension Band (2-piece)	_____	4-bolt Adapter Plate, 4 short, 4 long screws
_____	Neoplush Pelvic Band (Amputated side)	_____	Post-op Foot: L <input type="checkbox"/> R <input type="checkbox"/> Size: _____
_____	Neoplush Pelvic Band (Contralateral side)	_____	Titanium SACH Foot Adapter
_____	Per Socket Strap - All Sizes	_____	Titanium Rotatable 4-hole Connector
_____	Velcro® brand compatible neoprene socket band (all types)	_____	Titanium Pylon Tube
_____	Distal End Pad – per pkg. Size: _____	_____	S.S. SACH Foot Adapter
_____	Distal End Pad – Reticulated Size: _____	_____	S.S. Tube Clamp
_____	Aluminum Manual Locking Knee (<300 lbs.)	_____	S.S. Pylon Tube



**FLO-TECH® O & P Systems, Inc.**  
 PO Box 462 - 7325 Halseyville Road  
 Trumansburg, New York 14886  
 www.1800FLO-TECH.com  
 All materials © FLO-TECH® Modified March 2017

**Phone: 1-800-FLO-TECH** (356-8324)  
**Phone: (607) 387-3070**  
**Fax: (607) 387-3176**  
**E-mail: info@1800FLO-TECH.com**