

**Adult 7" Length Trans-Tibial Measurement Form**

Date: \_\_\_\_\_

PATIENT: \_\_\_\_\_ PROSTHETIST: \_\_\_\_\_ PO NUMBER: \_\_\_\_\_

**MEASUREMENTS** (ALL MEASUREMENTS SHOULD BE TAKEN OVER BANDAGES, SOCKS, ETC.): **RESIDUUM: Left \_\_\_\_\_ Right \_\_\_\_\_**  
**LENGTH MPT TO DISTAL END:** \_\_\_\_\_ **in/cm** (shorter than 4" or longer than 8¼" without bandages requires custom-made sockets)

**CIRCUMFERENCES: MPT \_\_\_\_\_ in/cm 2" UP FROM DISTAL END \_\_\_\_\_ in/cm** (Choose size based upon distal circ.)

*Sockets are automatically sized to provide the ideal anatomic knee circumference for the socket size ordered. Custom: \_\_\_\_\_*

Use the following sizing information for selecting APOPPS® sockets for your patient. Each of these sockets is a one-piece, polyethylene, clamshell design with anatomically trimmed overlapping sides, ideal for: accommodating bulbous, straight or conical shapes and gently reducing the residuum throughout the rehabilitation period and for weight bearing when used in conjunction with the UFOS™.

- ☞ The FLO-TECH-TOR™ is used as a prosthetic socket and for early supervised weight bearing prosthetic management of trans-tibial amputations.
- ☞ The VCSPS™ in conjunction with the UFOS™ is used for preparatory prosthetic management of trans-tibial amputations when unlimited knee motion is desired.
- ☞ The UFOS™ is an external frame outer socket used to attach the pylon and foot system. It is sized to match the FLO-TECH-TOR™ or VCSPS™ sockets. This allows your patient to wear the FLO-TECH-TOR™ or VCSPS™ sockets alone when the pylon and foot are not needed.

*The APOPPS® Program is designed to maintain a healthy amputation throughout the rehabilitation process.*

ADULT SIZES: 7" (4 ¼" to 6 ¼" LENGTH)	FLO-TECH-TOR™	VCSPS™	UFOS™	REHAB	PREP	APOPPS®
<b>#1416 – (MPT 14" to 16")</b> Distal: 13" (33cm) to 15" (38cm) Range: 11" (28cm) to 16" (41cm)	7-L _____ 7-R _____	7-L _____ 7-R _____	7-L _____ 7-R _____	7-L _____ 7-R _____	7-L _____ 7-R _____	7-L _____ 7-R _____
<b>#1618 – (MPT 16" to 18")</b> Distal: 15" (38cm) to 17" (43cm) Range: 13" (33cm) to 18" (46cm)	7-L _____ 7-R _____	7-L _____ 7-R _____	7-L _____ 7-R _____	7-L _____ 7-R _____	7-L _____ 7-R _____	7-L _____ 7-R _____
<b>#1821 – (MPT 18" to 21")</b> Distal: 17" (43cm) to 19" (48.5cm) Range: 15" (38cm) to 22" (56cm)	7-L _____ 7-R _____	7-L _____ 7-R _____	7-L _____ 7-R _____	7-L _____ 7-R _____	7-L _____ 7-R _____	7-L _____ 7-R _____
<b>#2124 – (MPT 21" to 24")</b> Distal: 20" (51cm) to 23" (58.5cm) Range: 17" (43cm) to 25" (63.5cm)	7-L _____ 7-R _____	7-L _____ 7-R _____	7-L _____ 7-R _____	7-L _____ 7-R _____	7-L _____ 7-R _____	7-L _____ 7-R _____

(9" length available with separate form. All other dimensions & lengths are custom orders)

**Foot & Pylon Units Available:** Includes pyramid plate, S.A.C.H. foot adapter, tube clamp and  
**Post-Op Wayfarer Foot** L R Size: 24 26 28 (Required with all foot orders)  
 Titanium Pylon Unit:  Stainless Steel Pylon Unit:

**OPTIONS & COMPONENT ITEMS:**

<u>QTY</u>	<u>DESCRIPTION</u>	<u>QTY</u>	<u>DESCRIPTION</u>
_____	Modification to Standard Socket	_____	4-bolt Adapter Plate, 4 short, 4 long screws
_____	Fork/Suspension Strap	_____	Titanium SACH Foot Adapter
_____	Per Socket Strap All Sizes	_____	Titanium Pyramid Socket Attachment Plate
_____	Velcro® brand Compatible neoprene socket band _____	_____	Titanium Pylon Tube w/Adapter
_____	Distal End Pad – per pkg. Size: _____	_____	Titanium Tube Clamp
_____	Distal End Pad Reticulated Size: _____	_____	Post-op Foot: <input type="checkbox"/> L <input type="checkbox"/> R Size: _____
_____	TT Waist belt/Suspension Belt	_____	Stainless Steel SACH Foot Adapter
_____	TT Padded/Lined Waist Belt	_____	Stainless Steel Pyramid Socket Attachment Plate
_____	Each Amputation Sock Size: _____	_____	Stainless Steel Pylon Tube
_____	Pour & Modify Cast	_____	Stainless Steel Tube Clamp