

RESIDUUM: LEFT ____ RIGHT ____ 29 DEGREES OR LESS ____ 30 DEGREES OR MORE ____

MEASUREMENTS (ALL MEASUREMENTS TAKEN OVER BANDAGES, SOCKS, ETC.):

LENGTH: MPT TO DISTAL END _____ In. (Cm) (shorter than 4" or longer than 8 1/4" without bandages requires custom made sockets)

CIRCUMFERENCES: MPT ____ In. (Cm) 2" UP FROM DISTAL END ____ In. (Cm) (Choose size based on this Circ.)
 Sockets are automatically sized to provide the ideal anatomic knee circumference for the socket size ordered.

ADULT SIZES

#1416 – (MPT 14" TO 16")
 Distal: 13" (33cm) to 15" (38cm)
 Range: 11" (28cm) to 16" (41cm)

MODIFIED FLO-TECH-TOR™

L _____ 7 or 9
 R _____

29 DEGREES OR LESS:

The anterior distal portion of a *FLO-TECH-TOR™* follows the angle of the contractures (up to 29 degrees) while the posterior section allows for graduated reduction of the contracture.

#1618 – (MPT 16" TO 18")
 Distal: 15" (38cm) to 17" (43cm)
 Range: 13" (33cm) to 18" (46cm)

L _____ 7 or 9
 R _____

#1821 – (MPT 18" TO 21")
 Distal: 17" (43cm) to 19" (48.5cm)
 Range: 15" (38cm) to 22" (56cm)

L _____ 7 or 9
 R _____

***30 DEGREES OR MORE:**

A modified thigh section of a *FLO-TECH-TOR™* & complete VCSPS™ with jointed uprights.

#2124 – (MPT 21" TO 24")
 Distal: 20" (51cm) to 23" (58.5cm)
 Range: 17" (43cm) to 25" (63.5cm)

L _____ 7 or 9
 R _____

These socket systems include:

- An additional neoprene band with Velcro® type attachment at the knee
- A long nylon stockinet pull-sock to allow for periodic re-positioning of tissue
- Slots cut in the distal pads
- Distal hole in the socket

*Create a jointed APOPPS® (Rehab and Prep prosthesis all in one) by adding an UFOS™ and waist belt, which allows for the attachment of Pylon and Foot components for ambulatory training. Remove the joints and thigh section (from the flexion control socket) when the contracture is eliminated to create a Prep System.

OPTIONS & COMPONENT ITEMS:

| <u>QTY</u> | <u>DESCRIPTION</u> | <u>QTY</u> | <u>DESCRIPTION</u> |
|------------|---|------------|---|
| _____ | Modification to Standard Socket | _____ | 4-bolt Adapter Plate, 4 short, 4 long screws |
| _____ | Fork/suspension strap | _____ | Titanium SACH Foot Adapter |
| _____ | Per Socket Strap - All Sizes | _____ | Titanium Pyramid Socket Attachment Plate |
| _____ | Velcro® Compatible Neoprene socket band _____ | _____ | Titanium Pylon Tube w/Adapter |
| _____ | Distal End Pad – per pkg. Size: _____ | _____ | Titanium Tube Clamp |
| _____ | Distal End Pad Reticulated Size: _____ | _____ | Post-op Foot: <input type="checkbox"/> L <input type="checkbox"/> R Size: _____ |
| _____ | TT Waist belt/suspension strap | _____ | Stainless Steel SACH Foot Adapter |
| _____ | TT Padded/Lined Waist belt | _____ | Stainless Steel Pyramid Socket Attachment Plate |
| _____ | Each Amputation Sock Size: _____ | _____ | Stainless Steel Pylon Tube w/Adapter |
| _____ | Pour & Modify Cast | _____ | Stainless Steel Tube Clamp |