

SOFTIE™ and SOFTIE™-TF Measurement Form

Date: _____

PATIENT: _____ PROSTHETIST: _____ PO NUMBER: _____

Trans-Tibial

1. A modified FLO-TECH-TOR™ made from 3/16" firm Pelite™ with Velcro®
2. An anatomically correct support strut, with neoprene attachment cover
3. A FLO-TECH® neoprene thigh band
4. A neoprene knee band
5. A FLO-TECH® neoprene socket band
6. A neoprene distal socket retainer (cap)
7. Standard 1" and 2" distal end pads

Regular

10" length (6¾" to 10")

Distal Range: 11" (28cm) to 23" (58.5cm)

L _____ R _____

Extra Large

10" length (6¾" to 10")

Distal Range: 20" (51cm) to 28" (71cm)

L _____ R _____

Trans-Femoral

1. A modified FLO-TECH-TOR™-TF made from 3/16" firm Pelite™ with Velcro®
2. A 3-piece pelvic band:
 - a. A 9" lateral socket extension band (used to extend the lateral trochanter area of the socket to interface with the trochanter section of the pelvic band)
 - b. A 2" by 5" double-sided Velcro® pressed between a standard neoprene trochanter section of the pelvic belt and the lateral socket extension band
 - c. The lateral socket extension band may be used either horizontally or vertically depending on the size of the patient.
 - d. It is attached to the SOFTIE™-TF with 2 patches of Velcro® at the proximal lateral area
3. A FLO-TECH® neoprene socket band
4. A neoprene distal socket retainer (cap)
5. 1 set 1" and 2" distal end pads

Regular:

13" length (9¾" to 13")

Distal Range: 13" (33cm) to 28" (71cm)

Universal Side QTY _____

Extra Large

13" length (9¾" to 13")

Distal Range: 23" (58.5cm) to 38" (96.5cm)

Universal Side QTY _____

OPTIONS AND REPLACEMENT ITEMS PRICES:

QTY	DESCRIPTION	TRANSTIBIAL
_____	Waist belt/Suspension strap	
_____	Padded & Lined Waist Belt/Suspension strap	
_____	SOFTIE modified fork strap	
_____	Neoprene Thigh Band	
_____	Neoprene Knee Band	
_____	Neoprene Distal Socket Band	
_____	Posterior Strut Neoprene Cover	
_____	Posterior Strut (only)	

ALL BANDS ARE VELCRO® COMPATABLE

QTY	DESCRIPTION	TRANSFEMORAL
_____	Neoprene Pelvic Suspension Band (3-Piece)	

_____	Neoprene Distal Retainer	
_____	Distal end pads/Per Pkg.	SIZE: _____
_____	Distal end pad/Reticulated	SIZE: _____
_____	Each Amputation Sock	SIZE: _____

FLO-TECH® O&P Systems, Inc.
PO Box 462 – 7325 Halseyville Road
Trumansburg, NY 14886

www.1800flo-tech.com

All materials © FLO-TECH® August 2017

Phone: 1-800-FLO-TECH (356-8324)

Phone: (607) 387-3070

Fax: (607) 387-3176

Email: info@1800flo-tech.com