

APOPPS®-TF

Trans-Femoral Measurement Form

Date: _____

PATIENT: _____ PROSTHETIST: _____ PO NUMBER: _____

Use the following sizing information for selecting APOPPS®-TF sockets for your patient. Each of these sockets is a one-piece, polyethylene, clamshell design with anatomically trimmed overlapping sides, ideal for accommodating bulbous, straight or conical shapes and gently reducing the residuum throughout the rehabilitation period.

- ⊗ The FLO-TECH-TOR™-TF is used as a prosthetic socket and for early supervised weight bearing prosthetic management of trans-femoral amputations.
- ⊗ The UFOS™-TF is an external frame outer socket used to attach the pylon and foot unit. It is sized to match the FLO-TECH-TOR™-TF you select for your patient.

The APOPPS® Program is designed to maintain a healthy amputation throughout the rehabilitation process.

MEASUREMENTS (ALL MEASUREMENTS SHOULD BE TAKEN OVER BANDAGES, SOCKS, ETC.): **RESIDUUM:** Left _____ Right _____

LENGTH Ischium to Distal End: _____ in/cm **Custom** _____

CIRCUMFERENCES:** Ischial Level _____ in/cm *3" Up from Distal End** _____ in/cm

**Top measurement at each size indicates ideal Ischial level circumference*

***Bottom measurement indicates the minimum to maximum circumference range at the Distal level.*

Sockets are automatically sized to fit the ideal circumference at the Above Ischial location.

12" (9¾" to 11¼" LENGTH)	FLO-TECH-TOR™	UFOS™	APOPPS®
SMALL – *16" (41cm) to 22" (56cm) **Range: 13" (33cm) to 17" (43cm)	L _____ R _____	L _____ R _____	L _____ R _____
MEDIUM – *19" (48cm) to 25" (63.5) **Range: 16" (41cm) to 20" (51cm)	L _____ R _____	L _____ R _____	L _____ R _____
LARGE – *22" (56cm) to 28" (71cm) **Range: 19" (48cm) to 23" (58.5cm)	L _____ R _____	L _____ R _____	L _____ R _____
X-LARGE – *26" (66cm) to 32" (81.5cm) **Range: 24" (61cm) to 28" (71.5cm)	L _____ R _____	L _____ R _____	L _____ R _____

Knee, Foot & Pylon Units Available: Includes manual locking knee, rotatable 4-hole connector, S.A.C.H. foot adapter, pylon and Post-Op Wayfarer Foot L R Size: 24 26 28 (Required with all foot orders)

Titanium/Knee/Pylon Unit Stainless Steel/Knee/Pylon Unit

OPTIONS & COMPONENT ITEMS:

QTY	DESCRIPTION	QTY	DESCRIPTION
_____	Modification to Standard Socket	_____	4-bolt Adapter Plate, 4 short, 4 long screws
_____	Neoprene distal socket band	_____	Titanium SACH Foot Adapter
_____	Per Socket Strap - All Sizes	_____	Titanium Rotatable 4-hole Connector
_____	Distal End Pad – per pkg. Size: _____	_____	Titanium Pylon Tube w/Adapter
_____	Distal End Pad Reticulated Size: _____	_____	Aluminum Manual Locking Knee (<300lbs)
_____	Complete Belt System	_____	Post-op Foot: <input type="checkbox"/> L <input type="checkbox"/> R Size: _____
_____	Each Amputation Sock Size: _____	_____	Stainless Steel SACH Foot Adapter
_____	Pour & Modify Cast	_____	Stainless Steel Rotatable 4-hole Connector
		_____	Stainless Steel Pylon Tube w/Adapter

FLO-TECH® O&P Systems, Inc.
PO Box 462 – 7325 Halseyville Road
Trumansburg, NY 14886
www.1800flo-tech.com
All materials © FLO-TECH® August 2017

Phone: 1-800-FLO-TECH (356-8324)
Phone: (607) 387-3070
Fax: (607) 387-3176
Email: info@1800flo-tech.com