

# APOPPS®-TF

## Trans-Femoral Measurement Form

Date: \_\_\_\_\_

PROSTHETIST: \_\_\_\_\_ PO NUMBER: \_\_\_\_\_

Use the following sizing information for selecting APOPPS®-TF sockets for your patient. Each of these sockets is a one-piece, polyethylene, clamshell design with anatomically trimmed overlapping sides, ideal for accommodating bulbous, straight or conical shapes and gently reducing the residuum throughout the rehabilitation period.

- ⊗ The FLO-TECH-TOR™-TF is used as a prosthetic socket and for early supervised weight bearing prosthetic management of trans-femoral amputations.
- ⊗ The UFOS™-TF is an external frame outer socket used to attach the pylon and foot unit. It is sized to match the FLO-TECH-TOR™-TF you select for your patient.

*The APOPPS® Program is designed to maintain a healthy amputation throughout the rehabilitation process.*

**MEASUREMENTS** (ALL MEASUREMENTS SHOULD BE TAKEN OVER BANDAGES, SOCKS, ETC.): **RESIDUUM:** Left \_\_\_\_\_ Right \_\_\_\_\_

**LENGTH Ischium to Distal End:** \_\_\_\_\_ in/cm **Custom** \_\_\_\_\_ \*

*\*Custom orders are not returnable or refundable, except for manufacturer defect.*

\*CIRCUMFERENCES: Ischial Level \_\_\_\_\_ in/cm \*\*3" Up from Distal End \_\_\_\_\_ in/cm

*\*Top measurement at each size indicates ideal Ischial level circumference*

*\*\*Bottom measurement indicates the minimum to maximum circumference range at the Distal level.*

*Sockets are automatically sized to fit the ideal circumference at the Above Ischial location.*

12" (9¾" to 11¼" LENGTH)	FLO-TECH-TOR™	UFOS™	APOPPS®
<b>SMALL – *16" (41cm) to 22" (56cm)</b> **Range: 13" (33cm) to 17" (43cm)	L _____ R _____	L _____ R _____	L _____ R _____
<b>MEDIUM – *19" (48cm) to 25" (63.5)</b> **Range: 16" (41cm) to 20" (51cm)	L _____ R _____	L _____ R _____	L _____ R _____
<b>LARGE – *22" (56cm) to 28" (71cm)</b> **Range: 19" (48cm) to 23" (58.5cm)	L _____ R _____	L _____ R _____	L _____ R _____
<b>X-LARGE – *26" (66cm) to 32" (81.5cm)</b> **Range: 24" (61cm) to 28" (71.5cm)	L _____ R _____	L _____ R _____	L _____ R _____

*Returns require authorization and are subject to restocking fee.*

### OPTIONS & COMPONENT ITEMS:

QTY	DESCRIPTION	QTY	DESCRIPTION
_____	Modification to Standard Socket	_____	4-bolt Adapter Plate, 4 short, 4 long screws
_____	Neoprene distal socket band	_____	Titanium SACH Foot Adapter
_____	Per Socket Strap - All Sizes	_____	Titanium Rotatable 4-hole Connector
_____	Distal End Pad – per pkg. Size: _____	_____	Titanium Pylon Tube w/Adapter
_____	Distal End Pad Reticulated Size: _____	_____	Aluminum Manual Locking Knee (<300lbs)
_____	Complete Belt System	_____	Post-op Foot: <input type="checkbox"/> L <input type="checkbox"/> R Size: _____
_____	Each Amputation Sock Size: _____	_____	Stainless Steel SACH Foot Adapter
_____	Pour & Modify Cast	_____	Stainless Steel Rotatable 4-hole Connector
		_____	Stainless Steel Pylon Tube w/Adapter

FLO-TECH® O&P Systems, Inc.  
PO Box 462 – 7325 Halseyville Road  
Trumansburg, NY 14886  
[www.1800flo-tech.com](http://www.1800flo-tech.com)  
All materials © FLO-TECH® March 20, 2018

Phone: 1-800-FLO-TECH (356-8324)  
Phone: (607) 387-3070  
Fax: (607) 387-3176  
Email: [info@1800flo-tech.com](mailto:info@1800flo-tech.com)