

# SOFTIE™ and SOFTIE™-TF Measurement Form

Date: \_\_\_\_\_

PROSTHETIST: \_\_\_\_\_ PO NUMBER: \_\_\_\_\_

## Trans-Tibial

1. A modified FLO-TECH-TOR™ made from 3/16" firm Pelite™ with Velcro®
2. An anatomically correct support strut, with neoprene attachment cover
3. A FLO-TECH® neoprene thigh band
4. A neoprene knee band
5. A FLO-TECH® neoprene socket band
6. A neoprene distal socket retainer (cap)
7. Standard 1" and 2" distal end pads

### Custom SOFTIE™ Measurements:

Length MPT to Distal End \_\_\_\_\_

### Circumferences:

MPT \_\_\_\_\_ 2" up from Distal End \_\_\_\_\_

## Regular

10" length (6¾" to 10")

Distal Range: 11" (28cm) to 23" (58.5cm)

L \_\_\_\_\_

R \_\_\_\_\_

## Extra Large

10" length (6¾" to 10")

Distal Range: 20" (51cm) to 28" (71cm)

L \_\_\_\_\_

R \_\_\_\_\_

## Trans-Femoral

1. A modified FLO-TECH-TOR™-TF made from 3/16" firm Pelite™ with Velcro®
2. A 9" lateral socket extension band (used to extend the lateral trochanter area of the socket to interface with the trochanter section of the pelvic band)
  - a. A 2" by 5" double-sided Velcro® pressed between a standard neoprene trochanter section of the pelvic belt and the lateral socket extension band
  - b. The lateral socket extension band may be used either horizontally or vertically depending on the size of the patient.
  - c. It is attached to the SOFTIE™-TF with 2 patches of Velcro® at the proximal lateral area
3. A FLO-TECH® neoprene socket band
4. A neoprene distal socket retainer (cap)
5. 1 set 1" and 2" distal end pads

## Regular:

13" length (9¾" to 13")

Distal Range: 13" (33cm) to 28" (71cm)

Universal Side

QTY \_\_\_\_\_

## Extra Large

13" length (9¾" to 13")

Distal Range: 23" (58.5cm) to 38" (96.5cm)

Universal Side

QTY \_\_\_\_\_

### OPTIONS AND REPLACEMENT ITEMS:

QTY	DESCRIPTION	TRANSTIBIAL
_____	TT Waist belt/Suspension strap	
_____	Padded & Lined Waist Belt/Suspension strap	
_____	SOFTIE modified fork strap	
_____	Neoprene Thigh Band	
_____	Neoprene Knee Band	
_____	Neoprene Distal Socket Band	
_____	Posterior Strut Neoprene Cover	
_____	Posterior Strut (only)	

### ALL BANDS ARE VELCRO® COMPATABLE

QTY	DESCRIPTION	TRANSFEMORAL
_____	Neoprene Pelvic Suspension Band (2-Piece)	
	(Required for TF)	

_____	Neoprene Distal Retainer	
_____	Distal end pads/Per Pkg.	SIZE: _____
_____	Distal end pad/Reticulated	SIZE: _____
_____	Each Amputation Sock	SIZE: _____

Returns require authorization and are subject to restocking fee.

FLO-TECH® O&P Systems, Inc.  
PO Box 462 – 7325 Halseyville Road  
Trumansburg, NY 14886

[www.1800flo-tech.com](http://www.1800flo-tech.com)

All materials © FLO-TECH® March 20, 2018

Phone: 1-800-FLO-TECH (356-8324)

Phone: (607) 387-3070

Fax: (607) 387-3176

Email: [info@1800flo-tech.com](mailto:info@1800flo-tech.com)