

PROSTHETIST: _____ PO NUMBER: _____

RESIDUUM: LEFT ____ **RIGHT** ____ **29 DEGREES OR LESS** ____ **30 DEGREES OR MORE** ____

MEASUREMENTS (ALL MEASUREMENTS TAKEN OVER BANDAGES, SOCKS, ETC.):

LENGTH: MPT TO DISTAL END ____ In. (Cm) (shorter than 4" or longer than 8 1/4" without bandages requires custom made sockets)

CIRCUMFERENCES: MPT ____ In. (Cm) **2" UP FROM DISTAL END** ____ In. (Cm) (Choose size based on this Circ.)

Sockets are automatically sized to provide the ideal anatomic knee circumference for the socket size ordered.

ADULT SIZES

#1416 – (MPT 14" TO 16")
Distal: 13" (33cm) to 15" (38cm)
 Range: 11" (28cm) to 16" (41cm)

MODIFIED FLO-TECH-TOR™

L ____ **7 or 9**
 R ____

29 DEGREES OR LESS:

The anterior distal portion of a *FLO-TECH-TOR™* follows the angle of the contractures (up to 29 degrees) while the posterior section allows for graduated reduction of the contracture.

#1618 – (MPT 16" TO 18")
Distal: 15" (38cm) to 17" (43cm)
 Range: 13" (33cm) to 18" (46cm)

L ____ **7 or 9**
 R ____

#1821 – (MPT 18" TO 21")
Distal: 17" (43cm) to 19" (48.5cm)
 Range: 15" (38cm) to 22" (56cm)

L ____ **7 or 9**
 R ____

***30 DEGREES OR MORE:**

A modified thigh section of a *FLO-TECH-TOR™* & complete VCSPS™ with jointed uprights.

#2124 – (MPT 21" TO 24")
Distal: 20" (51cm) to 23" (58.5cm)
 Range: 17" (43cm) to 25" (63.5cm)

L ____ **7 or 9**
 R ____

These socket systems include:

- An additional neoprene band with Velcro® type attachment at the knee
- A long nylon stockinet pull-sock to allow for periodic re-positioning of tissue
- Slots cut in the distal pads
- Distal hole in the socket

*Create a jointed APOPPS® (Rehab and Prep prosthesis all in one) by adding an UFOS™ and waist belt, which allows for the attachment of Pylon and Foot components for ambulatory training. Remove the joints and thigh section (from the flexion control socket) when the contracture is eliminated to create a Prep System.

OPTIONS & COMPONENT ITEMS:

Returns accepted for manufacturer defect only.

<u>QTY</u>	<u>DESCRIPTION</u>
_____	Modification to Standard Socket
_____	Fork/suspension strap
_____	Per Socket Strap - All Sizes
_____	Velcro® Compatible Neoprene socket band _____
_____	Distal End Pad – per pkg. Size: _____
_____	Distal End Pad Reticulated Size: _____
_____	TT Waist belt/suspension strap
_____	TT Padded/Lined Waist belt
_____	Each Amputation Sock Size: _____
_____	Pour & Modify Cast

<u>QTY</u>	<u>DESCRIPTION</u>
_____	4-bolt Adapter Plate, 4 short, 4 long screws
_____	Titanium SACH Foot Adapter
_____	Titanium Pyramid Socket Attachment Plate
_____	Titanium Pylon Tube w/Adapter
_____	Titanium Tube Clamp
_____	Post-op Foot: <input type="checkbox"/> L <input type="checkbox"/> R Size: _____
_____	Stainless Steel SACH Foot Adapter
_____	Stainless Steel Pyramid Socket Attachment Plate
_____	Stainless Steel Pylon Tube w/Adapter
_____	Stainless Steel Tube Clamp