

The FLO-TECH-TOR™ Provides:

*For peace of mind,
send your patients
home safely –
Send them home
with a
FLO-TECH-TOR™*



☞ **Minimizes the effects of potential injury during the early stages of rehabilitation**

Patients are sent home as soon as possible following an amputation. Many times patients either forget they have an amputation or are caught off guard while stumbling, causing the patient to land dangerously on the sutured or stapled wound. Even with heavy bandages, patients experience wound deprivation. These types of incidents are minimized through regular use of a FLO-TECH-TOR™.

☞ **Easy access to bandages and dressings for attentive wound care**

Sending patients home in a fiberglass or plaster cast makes wound care impossible, and leaves the patient in a heavy, awkward and uncomfortable device. The FLO-TECH-TOR™ is light and easy to remove and reapply, to allow the patient (or their attending aide) to provide diligent care to the amputation site.

☞ **Controlled early weight bearing, when the time is right**

Supervised controlled weight bearing can begin right away due to the flat bottom of the FLO-TECH-TOR™ providing increased stability. Controlled ambulation, under the supervision of a Physical Therapist, can begin without delay as well, when used with the UFOS™.

☞ **Control of swelling and edema**

There will come a time when the Physician wants some compression on the amputation. The FLO-TECH-TOR™ is ideal for this application and is removable to allow continued monitoring of elastic or non-elastic bandages and shrinker socks.

☞ **Maintenance of knee flexion and extension alignment**

Maintenance of good alignment of the knee is crucial to a successful rehabilitation program. The FLO-TECH-TOR™ is specifically designed to resist knee flexion contracture.

☞ **Volume control as it prepares the residual limb for a definitive prosthesis**

The FLO-TECH-TOR™ is fabricated of a semi-rigid plastic, which can be reduced in volume by adjusting the medial, lateral, anterior, posterior and circumference dimensions of the socket in the upper third, middle third or distal third of the socket as the healing progresses. This is accomplished by the socket's anatomic, overlapping design and one-piece construction; coupled with the use of a compressive outer band, use of shrinker socks or bandages inside the FLO-TECH-TOR™ to maintain total surface bearing and regular ambulation when used with the UFOS™.



FLO-TECH® O & P Systems, Inc.
PO Box 462 - 7325 Halseyville Road
Trumansburg, New York 14886
www.1800FLO-TECH.com
All materials © FLO-TECH® May 2019

Phone: 1-800-FLO-TECH (356-8324)
Phone: (607) 391-1070
Fax: (607) 391-1176
E-mail: info@1800FLO-TECH.com