

PROSTHETIST: _____ PO NUMBER: _____ Date: _____

MEASUREMENTS (ALL MEASUREMENTS SHOULD BE TAKEN OVER BANDAGES, SOCKS, ETC.): **RESIDUUM: Left _____ Right _____**

LENGTH MPT TO DISTAL END: _____ in/cm (shorter than 4" or longer than 8 1/4" without bandages requires custom-made sockets)

CIRCUMFERENCES: MPT _____ in/cm 2" UP FROM DISTAL END _____ in/cm (Choose size based upon distal circ.)

Sockets are automatically sized to provide the ideal anatomic knee circumference for the socket size ordered. **Custom:** _____

Use the following sizing information for selecting APOPPS® sockets for your patient. Each of these sockets is a one-piece, polyethylene, clamshell design with anatomically trimmed overlapping sides, ideal for: accommodating bulbous, straight or conical shapes and gently reducing the residuum throughout the rehabilitation period and for weight bearing when used in conjunction with the UFOS™.

- ☞ The FLO-TECH-TOR™ is used as a prosthetic socket and for early supervised weight bearing prosthetic management of trans-tibial amputations.
- ☞ The VCSPS™ in conjunction with the UFOS™ is used for preparatory prosthetic management of trans-tibial amputations when unlimited knee motion is desired.
- ☞ The UFOS™ is an external frame outer socket used to attach the pylon and foot system. It is sized to match the FLO-TECH-TOR™ or VCSPS™ sockets. This allows your patient to wear the FLO-TECH-TOR™ or VCSPS™ sockets alone when the pylon and foot are not needed.

The APOPPS® Program is designed to maintain a healthy amputation throughout the rehabilitation process.

ADULT SIZES: 7" (4 1/4" to 6 1/4" LENGTH)	FLO-TECH-TOR™	VCSPS™	UFOS™	REHAB	PREP	APOPPS®
#1416 – (MPT 14" to 16")	7-L _____	7-L _____	7-L _____	7-L _____	7-L _____	7-L _____
Distal: 13" (33cm) to 15" (38cm)	7-R _____	7-R _____	7-R _____	7-R _____	7-R _____	7-R _____
Range: 11" (28cm) to 16" (41cm)						
#1618 – (MPT 16" to 18")	7-L _____	7-L _____	7-L _____	7-L _____	7-L _____	7-L _____
Distal: 15" (38cm) to 17" (43cm)	7-R _____	7-R _____	7-R _____	7-R _____	7-R _____	7-R _____
Range: 13" (33cm) to 18" (46cm)						
#1821 – (MPT 18" to 21")	7-L _____	7-L _____	7-L _____	7-L _____	7-L _____	7-L _____
Distal: 17" (43cm) to 19" (48.5cm)	7-R _____	7-R _____	7-R _____	7-R _____	7-R _____	7-R _____
Range: 15" (38cm) to 22" (56cm)						
#2124 – (MPT 21" to 24")	7-L _____	7-L _____	7-L _____	7-L _____	7-L _____	7-L _____
Distal: 20" (51cm) to 23" (58.5cm)	7-R _____	7-R _____	7-R _____	7-R _____	7-R _____	7-R _____
Range: 17" (43cm) to 25" (63.5cm)						

(9" length available with separate form. All other dimensions & lengths are custom orders)

Foot & Pylon Units Available: Includes pyramid plate, S.A.C.H. foot adapter, tube clamp and Post-Op Wayfarer Foot L R Size: 24 26 28 (Required with all foot orders)

Titanium Pylon Unit Stainless Steel Pylon Unit

COMFORT OPTIONS & COMPONENT ITEMS:

QTY	DESCRIPTION	QTY	DESCRIPTION
_____	Modification to Standard Socket	_____	Each Amputation Socks Size: _____
_____	Waist belt/Suspension Belt	_____	Post-op Foot: L <input type="checkbox"/> R <input type="checkbox"/> Size: _____
_____	Fork/Suspension Strap	_____	Titanium SACH Foot Adapter
_____	Per Socket Strap - All Sizes	_____	Titanium Pyramid Socket Attachment Plate
_____	Velcro® brand Compatible neoprene socket band (all types)	_____	Titanium Pylon Tube w/Adaptor
_____	Distal End Pad – per pkg. Size: _____	_____	Titanium Tube Clamp
_____	Distal End Pad Reticulated Size: _____	_____	S.S. SACH Foot Adapter
_____	TT Waist belt/suspension strap	_____	S.S. Pyramid Socket Attachment Plate
_____	PADDED/LINED Waist Belt	_____	S.S. Pylon Tube w/Adaptor
_____	4-bolt Adapter Plate, 4 short, 4 long screws	_____	S.S. Tube Clamp



FLO-TECH® O & P Systems, Inc.
 PO Box 462 - 7325 Halseyville Road
 Trumansburg, New York 14886
 www.1800FLO-TECH.com
 All materials © FLO-TECH® Modified March 2017

Phone: 1-800-FLO-TECH (356-8324)
 Phone: (607) 387-3070
 Fax: (607) 387-3176
 E-mail: info@1800FLO-TECH.com